


SUBLIMITY ROUND SHAPE, SEMI-SMOOTH SURFACE



HIGH PROFILE - NATURGEL 1 - SOFT CONSISTENCY

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer*
85	68	35	52	LS 71 085	GAB 36 085
105	72	39	55	LS 71 105	GAB 36 105
125	78	40	58	LS 71 125	GAB 36 125
145	82	41	62	LS 71 145	GAB 36 145
165	86	43	65	LS 71 165	GAB 36 165
185	91	44	67	LS 71 185	GAB 36 185
205	93	45	70	LS 71 205	GAB 36 205
225	97	46	73	LS 71 225	GAB 36 225
245	99	47	75	LS 71 245	GAB 36 245
265	101	49	78	LS 71 265	GAB 36 265
285	103	50	80	LS 71 285	GAB 36 285
305	105	52	82	LS 71 305	GAB 36 305
335	109	54	85	LS 71 335	GAB 36 335
365	112	56	87	LS 71 365	GAB 36 365
395	114	57	89	LS 71 395	GAB 36 395
435	117	60	91	LS 71 435	GAB 36 435
475	122	61	95	LS 71 475	GAB 36 475
515	126	62	99	LS 71 515	GAB 36 515
565	129	63	102	LS 71 565	GAB 36 565
615	134	65	105	LS 71 615	GAB 36 615




*The sizers are prefilled with silicone gel with a smooth shell. They are delivered sterile, for single-use.

SUBLIMITY ROUND SHAPE, SEMI-SMOOTH SURFACE



SEMI-MODERATE PROFILE - NATURGEL 2 - FIRMER CONSISTENCY

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer*
80	82	21	56	LSC 76 080	GAB 30 080
100	87	23	59	LSC 76 100	GAB 30 100
120	92	24	62	LSC 76 120	GAB 30 120
140	95	26	65	LSC 76 140	GAB 30 140
160	103	27	69	LSC 76 160	GAB 30 160
180	107	28	71	LSC 76 180	GAB 30 180
200	111	29	73	LSC 76 200	GAB 30 200
220	114	30	76	LSC 76 220	GAB 30 220
245	118	31	78	LSC 76 245	GAB 30 245
270	122	32	80	LSC 76 270	GAB 30 270
295	126	33	83	LSC 76 295	GAB 30 295
325	130	34	85	LSC 76 325	GAB 30 325
370	136	35	88	LSC 76 370	GAB 30 370
415	142	36	91	LSC 76 415	GAB 30 415
470	148	38	94	LSC 76 470	GAB 30 470
510	153	39	98	LSC 76 510	GAB 30 510
560	157	42	102	LSC 76 560	GAB 30 560
-	-	-	-	-	-
-	-	-	-	-	-
-	-	-	-	-	-




The arc is the length from the highest point of projection to the implant base.



MODERATE PROFILE - NATURGEL 1 - SOFT CONSISTENCY

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer*
80	80	25	54	LS 90 080	GAB 35 080
100	84	28	59	LS 90 100	GAB 35 100
120	87	30	63	LS 90 120	GAB 35 120
140	91	31	66	LS 90 140	GAB 35 140
160	95	32	68	LS 90 160	GAB 35 160
180	101	33	72	LS 90 180	GAB 35 180
200	104	34	74	LS 90 200	GAB 35 200
220	108	35	75	LS 90 220	GAB 35 220
240	110	36	77	LS 90 240	GAB 35 240
260	112	38	78	LS 90 260	GAB 35 260
280	116	39	80	LS 90 280	GAB 35 280
300	119	40	82	LS 90 300	GAB 35 300
330	123	41	84	LS 90 330	GAB 35 330
360	126	42	87	LS 90 360	GAB 35 360
390	128	44	89	LS 90 390	GAB 35 390
430	133	45	93	LS 90 430	GAB 35 430
470	137	47	96	LS 90 470	GAB 35 470
510	140	49	98	LS 90 510	GAB 35 510
560	144	50	102	LS 90 560	GAB 35 560
610	147	51	108	LS 90 610	GAB 35 610




*The sizers are prefilled with silicone gel with a smooth shell. They are delivered sterile, for single-use.



HIGH PROFILE - NATURGEL 1 - SOFT CONSISTENCY

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer*
85	68	35	52	LS 91 085	GAB 36 085
105	72	39	55	LS 91 105	GAB 36 105
125	78	40	58	LS 91 125	GAB 36 125
145	82	41	62	LS 91 145	GAB 36 145
165	86	43	65	LS 91 165	GAB 36 165
185	91	44	67	LS 91 185	GAB 36 185
205	93	45	70	LS 91 205	GAB 36 205
225	97	46	73	LS 91 225	GAB 36 225
245	99	47	75	LS 91 245	GAB 36 245
265	101	49	78	LS 91 265	GAB 36 265
285	103	50	80	LS 91 285	GAB 36 285
305	105	52	82	LS 91 305	GAB 36 305
335	109	54	85	LS 91 335	GAB 36 335
365	112	56	87	LS 91 365	GAB 36 365
395	114	57	89	LS 91 395	GAB 36 395
435	117	60	91	LS 91 435	GAB 36 435
475	122	61	95	LS 91 475	GAB 36 475
515	126	62	99	LS 91 515	GAB 36 515
565	129	63	102	LS 91 565	GAB 36 565
615	134	65	105	LS 91 615	GAB 36 615




The arc is the length from the highest point of projection to the implant base.

SUBLIMITY ROUND SHAPE, SEMI-SMOOTH SURFACE

LS 74

SEMI-MODERATE PROFILE - NATURGEL 1 - SOFT CONSISTENCY

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer*
80	82	21	56	LS 74 080	GAB 33 080
100	87	23	59	LS 74 100	GAB 33 100
120	92	24	62	LS 74 120	GAB 33 120
140	95	26	65	LS 74 140	GAB 33 140
160	103	27	69	LS 74 160	GAB 33 160
180	107	28	71	LS 74 180	GAB 33 180
200	111	29	73	LS 74 200	GAB 33 200
220	114	30	76	LS 74 220	GAB 33 220
245	118	31	78	LS 74 245	GAB 33 245
270	122	32	80	LS 74 270	GAB 33 270
295	126	33	83	LS 74 295	GAB 33 295
325	130	34	85	LS 74 325	GAB 33 325
370	136	35	88	LS 74 370	GAB 33 370
415	142	36	91	LS 74 415	GAB 33 415
470	148	38	94	LS 74 470	GAB 33 470
510	153	39	98	LS 74 510	GAB 33 510
560	157	42	102	LS 74 560	GAB 33 560
-	-	-	-	-	-
-	-	-	-	-	-
-	-	-	-	-	-




*The sizers are prefilled with silicone gel with a smooth shell. They are delivered sterile, for single-use.

SUBLIMITY ROUND SHAPE, SEMI-SMOOTH SURFACE

LS 70

MODERATE PROFILE - NATURGEL 1 - SOFT CONSISTENCY

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer*
80	80	25	54	LS 70 080	GAB 35 080
100	84	28	59	LS 70 100	GAB 35 100
120	87	30	63	LS 70 120	GAB 35 120
140	91	31	66	LS 70 140	GAB 35 140
160	95	32	68	LS 70 160	GAB 35 160
180	101	33	72	LS 70 180	GAB 35 180
200	104	34	74	LS 70 200	GAB 35 200
220	108	35	75	LS 70 220	GAB 35 220
240	110	36	77	LS 70 240	GAB 35 240
260	112	38	78	LS 70 260	GAB 35 260
280	116	39	80	LS 70 280	GAB 35 280
300	119	40	82	LS 70 300	GAB 35 300
330	123	41	84	LS 70 330	GAB 35 330
360	126	42	87	LS 70 360	GAB 35 360
390	128	44	89	LS 70 390	GAB 35 390
430	133	45	93	LS 70 430	GAB 35 430
470	137	47	96	LS 70 470	GAB 35 470
510	140	49	98	LS 70 510	GAB 35 510
560	144	50	102	LS 70 560	GAB 35 560
610	147	51	108	LS 70 610	GAB 35 610



The arc is the length from the highest point of projection to the implant base.